



	QAS level - II
Application	Standardized, lyophilized monoclonal Anti-D IgG in specified concentration for Quality Assurance of reagents and techniques related to Matrix™ Gel System and classical tube technique
Reagent	Matrix™ QAS level -II contains: Standardized Anti-D IgG (Lyophilized)
Standardisation	Follows FDA guidelines
Storage / Stability	Lyophilized: 2-8°C, 12 months Reconstituted: 2-8°C, one week

Presentation	Pack	Cat.No.
Matrix™ QAS level -II	0.5 ml Anti-D and 10 ml 6% BSA	102780002



**MatrixTM
QAS Level-II****Intended Use:** _____

Standardized, lyophilized monoclonal Anti-D IgG in specified concentration for Quality Assurance of reagents and techniques related to MatrixTM Gel System and classical tube technique.

MatrixTM QAS Level-II kit components:

Reagent	0.5 ml Standardized Anti-D IgG (Lyophilized) and 10 ml 6% BSA
Other Accessories	Package Insert

FEATURES	BENEFITS
Standardized Anti-D IgG	Ideal for Quality Assurance of reagent and pipetting techniques
Lyophilized reagent	Highly stable

Storage / Stability	Temperature	Duration
Lyophilized	2-8°C	12 months
Reconstituted	2-8°C	One week

Available Pack Size
0.5 ml Anti-D and 10 ml 6% BSA