



ReaCell Panel

| | |
|----------------------------|---|
| Application | 0.8% reagent red blood cells for Antibody Screening on Matrix™ Gel Cards |
| Reagent | ReaCell Panel contains: <ul style="list-style-type: none"> • 11 cell panel of 0.8% reagent red blood cells for Antibody Screening on Matrix™ Gel Cards |
| Standardisation | Follows FDA guidelines |
| Storage / Stability | 2-8°C, 53 days |

| Presentation | Pack | Cat.No. |
|---------------|-----------|-----------|
| ReaCell Panel | 11 x 5 ml | IOTH-0241 |



ReaCell I II III

Antibody screening three cell panel for Matrix Gel Cards

Intended Use: _____

ReaCell I II III is 0.8% ready to use reagent red blood cells for Antibody Screening on Matrix Gel Cards.

ReaCell I II III Kit components:

| | |
|-------------------|---|
| ReaCell I II III | 3 cell panel each containing 0.8% red blood cells for Antibody Screening on Matrix Gel Cards. |
| Other Accessories | Antigen Chart, Package Insert |

| FEATURES | BENEFITS |
|---------------------|--|
| Convenient | Available in a hassle-free, easy to use and ready prepared solution Preservatives are in precise quantities ensure maximum shelf life of antigens |
| Single reagent pack | Having a shelf life of 53 days from the date of manufacturing and with added preservatives to maintain cell integrity |

| Storage / Stability | Temperature | Duration |
|---------------------|-------------|----------|
| Unopened kit | 2-8°C | 53 days |
| Opened kit | 2-8°C | 53 days |

| Available Pack Sizes |
|----------------------|
| 3 x 5 ml |