



	ANTI-A, ANTI-B, ANTI-A,B, ANTI-D
Application	Reagent combipack for ABO grouping, Rh ₀ (D) typing and confirmation by slide or tube test
Titre	<ul style="list-style-type: none"> • ANTI-A ≥ 1:256 • ANTI-B ≥ 1:256 • ANTI-A,B ≥ 1:256 • ANTI-D (IgM + IgG) ≥ 1:128
Class of Antibody	<ul style="list-style-type: none"> • ANTI-A, IgM, Murine Monoclonal • ANTI-B, IgM, Murine Monoclonal • ANTI-A,B, IgM, Murine Monoclonal • ANTI-D, IgM + IgG, Monoclonal Blend
Source	<ul style="list-style-type: none"> • ANTI-A, Murine hybridoma • ANTI-B, Murine hybridoma • ANTI-A,B, Murine hybridoma • ANTI-D, Human B cell line, EBV transformed
Specificity	<ul style="list-style-type: none"> • ANTI-A - 100% to A₁, A₂ and A_x antigens • ANTI-B - 100% to B antigens, Negative reacting with Acquired B characteristics • ANTI-A,B - 100% to A and B antigens, • ANTI-D - 100% to Rh₀(D) antigens
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
Reagent set for grouping typing and confirmation	4 x 10 ml	101950410
	4 x 12 ml	101950412





Eryscreen® Total

ABO / Rho(D) (IgM + IgG) reagent combipack for routine grouping and typing

Intended Use: _____

Eryscreen® Total, a reagent combipack for ABO grouping and Rho (D) (IgM+IgG) typing by slide or tube test.

Eryscreen® Total Kit components:

Monoclonal Total Anti A	Ready to use solution containing IgM (murine monoclonal) class antibodies specific to the 'A' antigen on the R.B.C
Monoclonal Total Anti B	Ready to use solution containing IgM (murine monoclonal) class antibodies specific to the 'B' antigen on the R.B.C
Monoclonal Total Anti D (IgM+IgG)	Ready to use solution containing IgM+IgG (human monoclonal) class antibodies specific to the 'Rho (D)' antigens on the R.B.C
Monoclonal Total Anti A,B	Ready to use solution containing IgM (murine monoclonal) class antibodies specific to the 'A' and 'B' antigens on the R.B.C
Other Accessories	Package Insert

Performance Characteristics:

Titre		Specificity: ANTI-A- 100% to A ₁ , A ₂ and A _x antigens ANTI-B - 100% to B antigens, negative reaction with Acquired B characteristics ANTI-D - 100% to Rho(D) antigens ANTI-A,B - 100% to 'A' and 'B' antigens, negative reaction with Acquired B characteristics Standardization: Follows AABB and FDA guidelines
Anti-A	≥ 1:256	
Anti-B	≥ 1:256	
Anti-A,B	≥ 1:256	
Anti-D (IgM+IgG)	≥ 1:256	

FEATURES	BENEFITS
Maintenance of optimum pH and ionic strength	Buffer with broad buffer capacity Ensures right conformational fit between antigen and antibody
Protect and preserve antibody activity	Optimum concentration of stabilizers Preserves antibody activity for long term use and storage
Optimized protein concentration	Well balanced to maintain optimal reactive conditions
Consistent reactivity on slide and tube tests	Reagent with high potential
Anti-A- 100% specific to A ₁ and A ₂ antigens Anti-B- 100% specific to B antigens Anti-D - 100% specific to Rho (D) antigens Anti-A,B - 100% specific to A and B antigens	Highly specific and truly negative with acquired B characteristics

Storage/Stability	Temperature	Duration
Unopened kit	2-8°C	24 Months
Opened kit	2-8°C	24 Months

Available Pack Sizes	
4 X 10 ml	4 X 12 ml