

R	EAGENS®	ReaCell I II III
Appl	ication	0.8% reagent red blood cells for Antibody Screening on Matrix™ Gel Cards
Reag	gent	ReaCell I II III contains:  • 3 cell panel of 0.8% reagent red blood cells

**Standardisation** 

Storage / Stability

Presentation	Pack	Cat.No.
ReaCell I II III	3 x 5 ml	IOTH-0240

2-8°C, 53 days

Follows FDA guidelines

for Antibody Screening on Matrix™ Gel Cards

## ReaCell I II III Antibody screening three cell panel for Matrix Gel Cards

Intended Use:	
intended USe:	

 $ReaCell\,I\,II\,III\,is\,0.8\%\,ready\,to\,use\,reagent\,red\,blood\,cells\,for\,Antibody\,Screening\,on\,Matrix\,Gel\,Cards.$ 

## ReaCell I II III Kit components:

ReaCell I II III	3 cell panel each containing 0.8% red blood cells for Antibody Screening on Matrix Gel Cards.
Other Accessories	Antigen Chart, Package Insert

FEATURES	BENEFITS
Convenient	Available in a hassle-free, easy to use and ready prepared solution Preservatives are in precise quantities ensure maximum shelf life of antigens
Single reagent pack	Having a shelf life of 53 days from the date of manufacturing and with added preservatives to maintain cell integrity

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	53 days
Opened kit	2-8°C	53 days

Available Pack Sizes	
3 x 5 ml	