

DH-615

Automatic Hematology Analyzer

Principles

Fluorescent nucleic acid staining for WBC differentiation
Impedance method for WBC, RBC and PLT test
Cyanide free colorimetry for HGB test

Throughput

Up to 100 tests per hour

Test Mode

CBC,CBC+RET,CBC+DIFF,CBC+DIFF+RET,RET

Parameters

35 reportable parameters + 29 research parameters

Sampling Mode

Auto sampling and closed sampling

Special Parameters

Immature granulocyte(IG) 、 Reticulocyte hemoglobin concentration (RET-HE)、
Immature platelet ratio (IPF%) 、 Immature reticulocyte ratio (IRF%)



Help Build Scalable
Smart Labs



10th Floor, Building B, High-tech Park, Guangqiao Road, Tianliao Community, Yutang Street, Guangming District, Shenzhen 518107, P. R. China

+86-755-26008015

www.dymind.com

+86-755-26008015

marketing@dymind.com



Declaration: Shenzhen Dymind Biotechnology Co., Ltd reserves the right to change the product of specifications and appearance at any time. For the information of this manual, Shenzhen Dymind Biotechnology Co., Ltd reserves the right to the interpretation and the decision.

P/N: EN-DH615[3.0]



DH-615

Automatic Hematology Analyzer



DH-615

Automatic Hematology Analyzer



Capillary blood sampling



Sample loading



Patient-reported outcomes



Recycling sample

The whole process of capillary blood measurement could be finished in 1-2 minutes.

Intelligent identification

- Supporting simultaneous batch measurement of whole blood and capillary blood.
- Automatically rotary to scan sample barcode, and manual adjustment is no need.

Intelligent measurement

- Automatically identify the type of sample and be able to mix capillary blood intelligently in every respect without breaking cells.
- Built-in capillary blood re-exam rules setting to help automatic re-examination.

Accurate Results

- Barcode management and standard mixing leads to more reliable results.
- Closed-cap measurement of capillary blood could avoid aerosol contamination and reduce the risk of cross-infection.

Leukocyte

providing immature granulocyte, atypical lymphocyte, high fluorescence intensity cells to reduce the missed detection rate of malignant hematological diseases;

Erythrocyte

providing LFR/MFR/HFR/RET-HE/IRF parameters to make precise diagnosis and monitor anemia diseases;



Two units can be connected to form a workstation to improve the inspection efficiency

Platelet

Dual-methodology PLT-O & PLT-I to avoid interferences from cell fragments, small red blood cells, platelet aggregation and etc.