


BRUNER

100 Strips

URS-4SG



REAGENT STRIPS FOR URINALYSIS

with Anti-Vc Ability

For in Vitro Diagnostic Use:

- Glucose
- Protein
- pH
- S • G

☐ URO

☐ LEU

☐ NIT

☒ PRO

☒ pH

☐ Ascorbat

☒ SG

☐ BLD

☐ KET

☐ BIL

☒ GLU

Important:

Store at temperature between 2-30°C(36-86 F) and out of direct sunlight. Replace cap immediately and tightly after taking out strips. Do not touch test areas. Do not desiccate from bottle. Read enclosed directions carefully.

Biosite
MEDICAL INSTRUMENTS

738 Shaw Boulevard, Mandaluyong
1554 Metro Manila

+63 9178512007 / +63 9178512008