



| ERYCARD™ ABO Blood Grouping Card Test for ABO Forward Grouping with Autocontrol | |
|---|---|
| Application | <ul style="list-style-type: none"> • Bed side ABO testing of Patient and Donor's blood unit before transfusion • Preconfirmation test, especially at nursing homes • For releasing compatible blood for emergency transfusions in <ul style="list-style-type: none"> • Defence Camps • Ambulance • Casualty services |
| Principle | Self performing, based on the principle of lateral flow guided by capillary action |
| Sensitivity & Specificity | 100% correlation with classical tube technique |
| Reagents | Reagents are pre-dried at the appropriate sample pad Anti-A (IgM) monoclonal antibodies in sample well A Anti-B (IgM) monoclonal antibodies in sample well B Autocontrol without any antibodies (Ctrl) Reagent Buffer (Proprietary Buffer) |
| Standardisation | Follows FDA guidelines |
| Test duration | Less than five minutes |
| Storage / Stability | 2-30°C, 24 months |

| Presentation | Pack | Cat.No. |
|--------------|------|---------|
| | | |

| | | |
|----------------------------------|----------|-----------|
| ERYCARD™ ABO kit contains | 24 Cards | 102840024 |
| Individually pouched device | | |
| Reagent buffer | | |
| Sample collection loops | | |





Erycard™ ABO

Blood Grouping Card for ABO Forward Grouping with Autocontrol

Intended Use: _____

Erycard™ ABO can be used as

- Bed side ABO & Rh testing of Patient and Donor's blood unit before transfusion
- Preconfirmation test, especially at nursing homes
- For releasing compatible blood for emergency transfusions in
 - Defence Camps
 - Ambulance
 - Casualty services

Erycard™ ABO Kit components:

| | |
|-------------------|--|
| Device | Pre-dried reagent at the appropriate sample pad Anti-A (IgM) monoclonal antibodies in sample well A Anti-B (IgM) monoclonal antibodies in sample well B Autocontrol without any antibodies (Ctrl) |
| Buffer | Ready to use Reagent Buffer (Proprietary Buffer) in a dropper bottle |
| Other Accessories | Desiccant pouch, Sample Collection Loops and Package Insert |

Performance Characteristics:

| | |
|---|---|
| Sensitivity & Specificity: 100% correlation with classical tube technique | Standardization: Follows FDA guidelines |
|---|---|

| FEATURES & BENEFITS |
|---------------------------------------|
| Blood Grouping in less than 5 minutes |
| Easy 3 step procedure |
| Simple interpretation of results |
| Auto control for test validation |
| Room temperature storage |

| Storage / Stability | Temperature | Duration |
|---------------------|-------------|-----------------|
| Unopened kit | 2-30°C | 24 months |
| Opened device | R.T. | Use immediately |

| Available Pack Sizes |
|----------------------|
| 24 Cards |