



	Diluent – 2 LISS
Application	For preparation of red blood cell suspensions used for Matrix™ gel system
Reagent	Buffered low ionic strength solution of appropriate sodium chloride molarity
Standardisation	Follows FDA guidelines
Storage / Stability	2-8°C, 18 months

Presentation	Pack	Cat.No.
Matrix™ Diluent-2 LISS	100 ml	102570100
	250 ml	10257250
	500 ml	102570500



**MatrixTM
Diluent-2 (LISS)**

Intended Use: _____

MatrixTM Diluent-2 (LISS) is Low ionic strength solution used for preparation of red cell suspension used for MatrixTM gel system.

MatrixTM Diluent-2 (LISS) kit components:

Matrix TM Diluent-2 (LISS)	Buffered low ionic strength solution of appropriate sodium chloride molarity.
Other Accessories	Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines

FEATURES	BENEFITS
Convenient	Ready to use Appropriate sodium chloride concentration to ensure that there is no hemolysis
Single Reagent Pack	With a shelf life of twelve months Is adequately buffered
Enhances Sensitivity	It enhances sensitivity by potentiating Ag-Ab reaction

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	12 Months
Opened bottle	2-8°C	12 Months

Available Pack Sizes	
100 ml	250 ml
500 ml	