## BRUNER LANCING DEVICE with Common Lancet

Lancing Device (lancet pen /pricker)

Easy to use

Durable

## **Product Description**

A device that holds a lancet firmly, and when triggered, moves the lancet linearly ahead to prick the skin in a controlled manner. Lancing devices usually can be adjusted to alter the depth that the lancet goes into the skin. The purpose is to obtain a small amount of capillary blood for testing.





Manila Office: 738 Shaw Boulevard, Greenfield District,

Mandaluyong City, Metro Manila

Tel: (02)7505-7425 / (02)7505-7443 Phone: (0917)-851-2007 / (0917)-851-2008

Davao Office: 553-555 Manga St., Juna Subdivision Matina, Davao City

www.biositemedicalinstruments.com.ph