



	ABO Subgrouping Card Anti-H
Application	For detection of H antigen
Card Configuration	6 x Anti-H / Six profiles per card
Reagent	Matrix™ ABO Subgrouping Card Anti-H contains: <ul style="list-style-type: none"> • Anti-H lectin
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ ABO Subgrouping Card Anti-H	24 Tests	102760024



**MatrixTM**
ABO Subgrouping Card Anti-H**Intended Use:** _____MatrixTM ABO Subgrouping Card Anti-H is used for detection of H antigens on RBCs.**MatrixTM ABO Subgrouping Card Anti-H kit components:**

MatrixTM ABO Subgrouping Card Anti-H	Matrix TM ABO Subgrouping Card Anti-H contains Anti-H lectin
Configuration	6 x Anti-H/ Six profiles per card
Other Accessories	Card fitment and Package Insert

Performance Characteristics:**Standardization:** Follows FDA guidelines

FEATURES	BENEFITS
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards