

Multi-type Sample DNA/RNA Extraction-Purification Kit (Magnetic Beads Method)



Application:

The magnetic beads and buffer system with unique separation effect can be used to extract high-purity viral DNA/RNA from samples quickly, highly sensitively and efficiently.

The extracted and purified nucleic acid can be used in various common downstream experiments such as restriction digestion, reverse transcription, PCR, RT-PCR, Southernblot, etc.

Method Advantage:

- ①. Safe and efficient
- ②. Less time requirement

Extraction process:

Add lyse- add magnetic bead- separate magnetic bead-washing- Separate magnetic

Product advantage:


- ①. High sensitivity: High-quality nucleic acid can be extracted after the positive samples are diluted by 104 and 107 respectively 10^4
- ②. Good repeatability: Good repeatability of extracting the same sample



Reagent For BNP32/48



Reagent For BK-HS32/AutoHS96



Reagent For BNP96



Reagent For BK-HS96

Parameters:

Name	Specification	Note Mark
Nucleic Acid Extraction Kit I	50T/Box	Manual
Nucleic Acid Extraction Kit II	32T/64T/96T/Box	BK-HS32/ BK-AutoHS96
Nucleic Acid Extraction Kit III	96T/Box	BK-HS96
Nucleic Acid Extraction Kit IV	96T/Box	BNP96
Nucleic Acid Extraction Kit V	64T/Box	BNP32&48