





SARS-CoV-2 RT-PCR Test Kit



## SARS-CoV-2 RT-PCR Test Kit

## Oualitative Detection of Nucleic Acid from SARS-CoV-2



CE marked, UK MHRA approved. Able to detect the UK variant (B.1.1.7), South Africa variant (B.1.351) and Brazil variant (P.1) mutant strains.



Triple-target assay in single well, covering highly conserved region within N, E and ORF1ab (RdRp) genes for highly specificity



811 clinical sample were used to verify the superiority



Full-process controls with negative, positive and internal controls, programmed time less than 60 min



Applied PCR machine:
Applied Biosystems® 7500, Applied
Biosystems® QuantStudio5, Promotor®
960 Real-Time PCR System, SLAN®-96P,
CFX96 (Bio-Rad), LightCycle® 480 (Roche)
Rotor-Gene Q(Qiagen), MA6000,
AGS9600, X960, Gentier 96,
Quant Gene9600

## **Key Features and Performance**

Technology	Real-Time reverse transcription polymerase chain reaction (rRT-PCR)		
Target Genes	N, E, ORF1ab (RdRp)		
Specimen Types	Oropharyngeal swab, nasopharyngeal swab and deep cough sputum		
Internal Control	RNase P		
Minimum Sample Volume	200 μL		
Programmed Duration	<60 minutes (Test duration depends on applied PCR machine)		
Limit of Detection	500 copies/mL		
Precision (In-batch and batch)	CV<5%		
Sensitivity	99.79% (95% CI: 98.82%-99.99%)		
Specificity	99.04% (95% CI: 98.11%-99.58%)		
Accuracy	99.31% (95% CI: 98.69%-99.68%)		
Storage Temperature	- 20 ± 5 °C		
Package Size	32 Tests/Kit		

## **Ordering Information**

Product Name	Catalog No.	Components	
SARS-CoV-2 RT-PCR Test Kit		Reaction Mix	1 x 630 μL
	P131-1591√	Enzyme Mix	1 x 70 μL
	F 131-1391V	Positive Control	1 x 200 μL
		Negative Control	1 x 200 μL

√CE marked

Programs and specifications are subject to change without prior notice. Please contact local distributors.



ACON Biotech (Hangzhou) Co., Ltd.

No.210 Zhenzhong Road, West Lake District, Hangzhou, P.R.China, 310030 Tel: +86-571-8777 5689 Fax: +86-571-8777 5781 Email: PCR.Info@aconlab.com.cn http://www.aconbio.com