

Diagnostic Kit for Coxsackievirus A16 RNA

(PCR-Fluorescence Probing)

CA16



Accurate Detection



EXTREMELY SIMPLE WORKFLOW

1. Upload sample and reagent

2. Select program

3. Detect and analyze



FEATURES







IC MONITOR

Lentivirus monitors the entire process

FULLY AUTOMATED

Matching Natch CS2, Extraction standardization

HIGH PRECISION

Advanced super-paramagnetic nano-beads specifically binding with sample RNA

PARAMETERS



Product Features	Parameter
Specimen Types	Throat swab
Technical Platform	Advanced magnetic beads technology
Target Genes	CA16 RNA
Internal Control	Lentivirus diluent
PCR Instruments	ABI7500 and Stratagene Mx3000P
Amplification Time	100min
Sensitivity	200 copies/mL

ORDER INFORMATION



No.	Product Name	Spec.
S3020E	Diagnostic Kit For Coxsackievirus A16 RNA (PCR-Fluorescence Probing)	24T/Kit



ADDR.: No. 680, Lusong Road, Yuelu District, Changsha, Hunan Province, 410205, P.R.China

TEL.: +86-731-88883176-6116 **FAX**: +86-731-89736985

MAIL: info@sansure.com.cn / support@sansure.com.cn

WEB: www.sansureglobal.com

