



ERYGEN-ABO	Red Blood Cell Panel for reverse grouping
Application	3% Reagent red blood cell panel for confirmation of the results of forward grouping
Specificity	ABO isoagglutinins
Reagent	Erygen-ABO panel of 3% reagent red blood cells comprises of: <ul style="list-style-type: none"> • Red blood cells of group A1 • Red blood cells of group B • Red blood cells of group O
Standardisation	Follows FDA guidelines
Storage / Stability	2-8°C, 50 days

Presentation	Pack	Cat.No.
Erygen – ABO	3x2ml	101950032





Erygen-ABO

Red Blood Cell panel for reverse grouping

Intended Use: _____

Erygen-ABO comprises of a 3% reagent red blood cells of A1, B and O suspended in an isotonic medium containing a preservative which is employed in the confirmation of the results of forward grouping.

Erygen-ABO Kit components:

Erygen-ABO - A ₁ cells	3% reagent red blood cells of Group A1
Erygen-ABO - B cells	3% reagent red blood cells of Group B
Erygen-ABO - O cells	3% reagent red blood cells of Group O
Other Accessories	Package Insert

Performance Characteristics:

Specificity: ABO isoagglutinins	Standardization: Follows FDA guidelines
--	--

FEATURES	BENEFITS
Ready to use	Solution comprises of RBCs in an appropriate concentration in a preservative solution Components in precise quantities ensure maximum benefits during blood grouping
Four component pack	Separate reagents for different cell types

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	50 days
Opened kit	2-8°C	50 days

Available Pack Sizes
3 X 2 ml