

Data sheet - Ice Spray

CE IVD	<i>In vitro diagnostic medical device</i>	Supplier <i>Histo-Line Laboratories Srl</i>
--------	---	---

<i>Code Product</i>	<i>Packaging</i>
R0020	12 bottles of 200 ml each
Code CND: W01030799	

Applications

Freezing spray for histology.

The product allows to quickly freeze tissues to be cut in the cryostat and to cool the paraffin embedded before microtome cutting where a cooling plate is not available. Spray directly on the sample to be cooled. Does not contain fluorinated greenhouse gases, so the Regulation (EU) 842/06 article 7 par 2 is not applicable.

Storage

Store in a cool and dry place, at room temperature and away from heat sources.

Do not expose the propellant gas to open flames.

Validity

Product validity: 5 years. It is advisable to always check the validity of the product on the packaging.

N.B. For anything not mentioned in this technical data sheet: compositions, precautionary measures, risk phrases, storage consult, fire-fighting measures and what else, consult the Safety Data Sheet downloadable from our website: www.histoline.com or request it at: specialist@histoline.com

Milan January 2018