

Technical Data Sheet PAP Pen

CE IVD	For In Vitro Diagnostic Use	Manufacturer: Histo-Line Laboratories Srl
--------	-----------------------------	---

Code Product	Packaging
LP 0001 Super PAP Pen	1 pc.
LP 0002 Mini PAP Pen	1 pc.



Intended Use

PAP Pen

This marking pen has been designed to provide a thin film-like barrier when a circle is drawn around the specimen on a slide. This barrier creates the proper surface tension to hold an antibody solution within the target area on the slide.

General features

The barrier creates the proper surface tension to hold an antibody solution within the target area on the slide. PAP Pen contains a special formulation which is insoluble in alcohol and acetone.

It can be removed, if desired, by xylene after the staining procedure is completed.

Validity

Validity is not applicable to this type of product.

Storage Condition

Store with cap well closed, in a fresh area.

N.B. Although not reported in this technical sheet: compositions, precautionary measures, risk phrases, refer to storage, fire precautions, and more, refer to the Safety Data Sheet can be downloaded from our website: www.histoline.com or request it to: info@histoline.com

Milan Jul 2022