



# BRUNER

100 Strips

# URS-10A



## REAGENT STRIPS FOR URINALYSIS

with Anti-Vc Ability

*For in Vitro Diagnostic Use:*

- Leukocytes
- Nitrite
- Urobilinogen
- Protein
- pH
- Blood
- Specific Gravity
- Ketone
- Bilirubin
- Glucose

## Important:

Store at temperature between 2-30°C(36-86 F) and out of direct sunlight. Replace cap immediately and tightly after taking out strips. Do not touch test areas. Do not desiccant from bottle. Read enclosed directions carefully.

**Biosite**  
MEDICAL INSTRUMENTS

738 Shaw Boulevard, Mandaluyong  
1554 Metro Manila

+63 9178512007 / +63 9178512008