



	<b>ABO Subgrouping Card Anti-A1</b>
<b>Application</b>	For confirmation of A1 subgroup of A
<b>Card Configuration</b>	6 x Anti-A1 / Six profiles per card
<b>Reagent</b>	Matrix™ ABO Subgrouping Card Anti-A1 contains: <ul style="list-style-type: none"> <li>• Anti-A1 lectin</li> </ul>
<b>Standardisation</b>	Follows FDA guidelines
<b>Storage / Stability</b>	4-25°C, 12 months

Presentation	Pack	Cat.No.
<b>Matrix™ ABO Subgrouping Card Anti-A1</b>	24 Tests	102750024





## Matrix™ ABO Subgrouping Card Anti-A<sub>1</sub>

**Intended Use:** \_\_\_\_\_

Matrix™ ABO Subgrouping Card Anti-A<sub>1</sub> is used for confirmation of A<sub>1</sub> subgroup of A.

### Matrix™ ABO Subgrouping Card Anti-A<sub>1</sub> kit components:

<b>Matrix™ ABO Subgrouping Card Anti-A<sub>1</sub></b>	Matrix™ ABO Subgrouping Card Anti-A <sub>1</sub> contains Anti-A <sub>1</sub> lectin
<b>Configuration</b>	6 x Anti-A <sub>1</sub> / Six profiles per card
<b>Other Accessories</b>	Card fitment and Package Insert

### Performance Characteristics:

**Standardization:** Follows FDA guidelines

FEATURES	BENEFITS
<b>Reliable and consistent reactivity only with A<sub>1</sub> and A<sub>1</sub>B cells</b>	Optimized to exclude A <sub>2</sub> reactivity and weaker subgroups Calibrated to pick up only A <sub>1</sub> and A <sub>1</sub> B cells
<b>Stable End Reactions</b>	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
<b>Objective Interpretation of results</b>	Easy Interpretation of weak reactions Relatively skill free
<b>Small sample size</b>	More flexible More convenient
<b>Indigenously Manufactured</b>	Easy availability Long shelf life Another first from Tulip
<b>Tulip's Technical Support</b>	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards