



	ERYGEN-ABO
Application	0.8% Reagent red blood cell panel for confirmation of the results of forward grouping
Specificity	ABO isoagglutinins
Reagent	Matrix Erygen-ABO panel of 0.8% reagent red blood cells comprises of: <ul style="list-style-type: none"> • Red blood cells of group A1 • Red blood cells of group B • Red blood cells of group O
Standardisation	Follows FDA guidelines
Storage / Stability	2-8°C, 36 days

Presentation	Pack	Cat.No.
Matrix™ ERYGEN-ABO	3 x 2 ml	102070032





Matrix™ Erygen-ABO Red Blood Cell panel for reverse grouping

Intended Use: _____

Matrix™ Erygen-ABO comprises of a 0.8% reagent red blood cells of A1, B and O suspended in an isotonic medium containing a preservative which is employed in the confirmation of the results of forward grouping.

Matrix™ Erygen-ABO Kit components:

Matrix™ Erygen-ABO - A ₁ cells	0.8% reagent red blood cells of Group A1
Matrix™ Erygen-ABO - B cells	0.8% reagent red blood cells of Group B
Matrix™ Erygen-ABO - O cells	0.8% reagent red blood cells of Group O
Other Accessories	Package Insert

Performance Characteristics:

Specificity: ABO isoagglutinins	Standardization: Follows FDA guidelines
--	--

FEATURES	BENEFITS
Ready to use	Solution comprises of RBCs in an appropriate concentration in a preservative solution Components in precise quantities ensure maximum benefits during blood grouping
Four component pack	Separate reagents for different cell types

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	36 days
Opened kit	2-8°C	36 days

Available Pack Sizes
3 X 2 ml