

S10 S12

Patient Monitor



Size and Weight

Size	S12:	175mm X 320mm X 262mm
	S10:	168mm X 288mm X 236mm
Weight	< 4kg	

Power

Standard According to IEC 60601-1 and IEC 60601-1-2

Input voltage	AC (100-240) V(±10%)
Frequency	50Hz/60Hz
Input power	100VA

Display

Type	Color TFT LCD
Size(diagonal)	12.1" / 10.4" (S12 / S10)
Resolution	800×600 pixels

Recorder(Optional)

Type	Thermal dot array
Paper width	50 mm ±1mm
Recording speed	12.5 mm/s, 25 mm/s, 50 mm/s
Recording waveform	Maximum 3 tracks

Battery

Type	Rechargeable Li-ion battery 11.1V 2.5Ah / 5.0Ah
Operating time	>240 / 480 minutes (2.5Ah / 5.0Ah)
(1 new and fully charged battery at 25°C temperature, connecting SpO2 sensor & NIBP work on AUTO mode for 30 minutes interval)	
Charge time	<6 / 12 hours(2.5Ah / 5.0Ah)

Data Storage

Alarm event	3000 groups and associated waveform
Trend	180h, minimum resolution is 1min
	6h, minimum resolution is 5s
ARR event	3000 groups and associated waveform
NIBP	2400 groups
Holographic waveform	72 hours

Interfacing & I/O devices

Keyboard & Mouse	Support
Barcode Scanner	Support 1D barcode (USB connector)
Wired network	1 standard RJ45 interfaces
Wifi (option)	Protocol: IEEE802.11a/b/g/n
Wifi frequency	Dual Band: 2.4G/5G
USB socket	2 sockets
Video output	1 VGA (option)
Multifunctional port	nurse call / defibrillation sync. / analog output

ECG

Lead	3 lead: I, II, III
	5 lead: I, II, III, aVR, aVL, aVF, Vx
	6-lead: I, II, III, aVR, aVL, aVF, Va, Vb
	12-lead: I, II, III, aVR, aVL, aVF, V1~V6 (S12 option)
	Auto: identify leads automatically

CMRR	Monitor / Operation mode ≥ 110 dB
	Diagnostic mode ≥ 100 dB
Bandwidth (-3dB)	Monitor mode: 0.5 Hz to 40 Hz
	Operation mode: 1 Hz to 25Hz
Input impedance	≥ 5.0 MΩ
Input signal range	-10.0mV~+10.0mV
Electrode offset potential	± 500 Mv d.c.
System noise	≤ 30 μVpp (RTI)
Recovery time after defibrillation:	waveform recover to baseline in 10s
Sweep speed	6.25mm/s, 12.5 mm/s, 25 mm/s, 50mm/s.

ST segment

Measurement range	-2.0 mV to +2.0 mV
Accuracy	-0.8 mV to +0.8 mV: ±0.02 mV or ±10% (whichever is greater)
Resolution	0.01mV

Heart Rate

Measurement range	Adult	10 bpm to 300 bpm
	Pediatric & Neonatal	10 bpm to 350 bpm
Resolution	1 bpm	
Accuracy	±1% or ±1 bpm, whichever is greater	

Arrhythmia analysis

27 Kinds (ASYSTOLE, BRADYCARDIA, TACHYCARDIA, EXTREME BRADYCARDIA, EXTREME TACHYCARDIA, VENTRICULAR BRADYCARDIA, VENTRICULAR TACHYCARDIA, NONSUSTAINED VENTRICULAR TACHYCARDIA, VENTRICULAR FIBRILLATION, ATRIAL FIBRILLATION, ATRIAL FIBRILLATION END, R ON T, VENTRICULAR RHYTHM, PNC, PNP, PAUSE, PVC, PAUSES/MIN HIGH, RUNNING PVCS, COUPLET, BIGEMINY, TRIGEMINY, FREQUENT PVCS, MISSED BEAT, ECG NOISE, IRREGULAR RHYTHM, IRREGULAR RHYTHM END)

Respiration

Lead	Selected from: I (RA-LA) or II (RA-LL)
Measurement range	0 rpm to 150 rpm
Resolution	1 rpm
Accuracy	±2 rpm or ±2% , whichever is the greater
Delay of apnea alarm	Adjustable delay time: 10s ~ 60s

QT analysis

Measurement range	QT: 200ms~700ms
	QTc: 200ms~700ms
	ΔQTc: -500ms~500ms
	QT-HR: Adult: 15bpm~150bpm
	Pediatric/neonatal:15bpm~180bpm
Resolution	QT, QTc, ΔQTc: 1ms
	QT-HR: 1bpm
Accuracy	QT: ±30ms

NIBP

Measurement way	Automatic oscillometry	
Measurement mode	Manual , Auto, STAT	
Intervals for Auto measurement: 1/2/2.5/3/5/10/15/20/30min, 1/1.5/2/4/8h		
STAT mode cycle time 5 minutes.		
Systolic range	Adult	30 to 270 mmHg
	Pediatric	30 to 235 mmHg
	Neonatal	30 to 135 mmHg
Diastolic range	Adult	10 to 220 mmHg
	Pediatric	10 to 220 mmHg
	Neonatal	10 to 110 mmHg
Mean range	Adult	20 to 235 mmHg

Pressure accuracy	Pediatric	20 to 235 mmHg
	Neonatal	20 to 125 mmHg
	Static:	±3 mmHg
	Clinic:	mean error ±5 mmHg
	Standard deviation:	≤8 mmHg
Inflation time for cuff	Less than 40s. (standard adult cuff)	
Cuff pressure range	0 to 300 mmHg	
PR range	40 bpm to 240 bpm	
Measurement time	20s to 45s (typical value)	
Lead standard	AHA, IEC	
Gain	Auto, 2.5 mm/Mv (×0.25), 5 mm/mV (×0.5), 10 mm/mV (×1), 20 mm/mV (×2), 40 mm/mV (×4)	

BLT SpO2

Measurement range	0% ~ 100%
Accuracy(clinical)	70% ~ 100% ≤3% (SpO2 probe included) 0% ~ 69% unspecified

PR

Measurement range	25 bpm to 300 bpm
Resolution	1bpm
Accuracy	± 3bpm

PI

Measurement range	0.05~20.00%
Resolution	0.01%
Accuracy	±0.1% or ±10% of reading, whichever is greater

RESP (from pleth)

Measurement range	0 rpm ~90 rpm
Resolution	1 rpm
Accuracy	± 2rpm

Temperature (Dual-Temp for S12 only)

Parameter	T1,T2,TD
Probe	YSI400 series probe (2252 Ω @25℃)
Measurement range	0.0℃ to 50.0℃(32℉ to 122℉)
Accuracy	±0.1℃ or ±1℉ (exclusive of probe)
Resolution	0.1℃ or 1℉
Unit	℃ or ℉

IBP (option for S12 only)

Sensitivity of transducer	5uV/V/ mmHg, ±2%
Impedance of transducer	300Ω to 3000Ω
Measurement range	-50 mmHg to +360 mmHg
Measurement accuracy	±2 mmHg or ±2% of the reading, whichever is the greater (exclusive of transducer)
Resolution	1 mmHg
Unit	mmHg, kPa, cmH2O
Transducer sites	ART/CVP/ICP/PA/Ao/UAP/BAP/FAP//LAP/RAP/UVPLV/PAWP, additionally, P1 & P2 are arbitrary sites

PPV

Measurement range	0~50%
Resolution	1.00%

PR

Measurement range	30 bpm to 300 bpm
Resolution	1bpm
Accuracy	±1% or ±1bpm whichever is greater

Software overpressure protection	Adult	(297±3) mmHg
	Pediatric	(252±3) mmHg
	Neonatal	(147±3) mmHg

MicroFlow CO2 (option for S12 only)

Measurement range	0% to 25% (0 mmHg to 190 mmHg)
Unit	0.1% or 1mmHg
Unit	%, mmHg, kPa
Accuracy	± (0.43% + 8% of reading)
Preheating time	<10s (Report concentration and achieve highest accuracy)
Rise time	<3s (including delay time and rise time)
Sample Flow Rate	50±10mL/min
awRR range	0 rpm to 150 rpm
awRR accuracy	±1 rpm

Mainstream CO2 (option for S12 only)

Measurement range	0% to 25% (0 mmHg to 190 mmHg)
Resolution	0.1% or 1mmHg
Preheating time	<10s
Rise time	<90ms
Unit	%, mmHg, kPa
Accuracy	± (0.43% + 8% of reading)
awRR range	0 rpm to 150 rpm
awRR accuracy	±1 rpm

C.O. (option for S12 only)

Measurement range	C.O.	0.1 L/min to 20 L/min
	TB	23.00℃ ~ 43.00℃
	TI	-1.0℃ ~ 27.0℃
Resolution	C.O.	0.1 L/min
	TB	0.01℃
	TI	0.1℃
Accuracy	C.O.	±5% or ±0.1L/min, whichever is greater
	TB	±0.1℃
	TI	±0.1℃

Drip Monitor (option)

Measurement range	Drip rate	5~200 Drops/min (1mL of conventional tube =20 drops)
	Accuracy	±2 digit or ±2% (whichever is greater)
Unit	Drops/min, mL/h, can be automatically converted (1mL conventional tube=20 drops is mainly used.)	
	Liquid stop function	Alarm and stop liquid when infusion is completed. Alarm when drip rate is abnormal.

Standard configuration:

3/5/6 lead ECG, HR, SpO2, PI, RESP(from pleth), NIBP, Temp, Dual-Temp(S12), Capacitive Touch Screen, Rechargeable Li-ion battery (2.5Ah).

Option:

S10: Drip monitor(DM), Rechargeable Li-ion battery (5Ah),

S12: Drip monitor(DM), 12 lead ECG, Voice assistant, Nurse call / Defibrillation sync. /VGA output, Rechargeable Li-ion battery (5Ah). 2-IBP, C.O., Mainstream/Microflow EtCO2.

Others: Thermal Printer, Rolling stand, Wall mount



*Specifications subject to change without prior notice.

Add: No.2 Innovation First Road, Technical Innovation Coast, Hi-tech Zone, Zhuhai, P.R.China
Tel: +86-756-3399935 3399999 Fax: +86-756-3399911
E-Mail: Overseas@blt.com.cn Post code: 519085
Http://www.blt.com.cn

S10/S12-210325-eng-V1; Date: Mar. 2021

