



 Forward and Reverse Grouping Card with Sub Grouping	<b>Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping</b>
<b>Application</b>	For determination of ABO and Rho (D) blood group systems with subgrouping and reverse grouping for confirmation of ABO grouping
<b>Card Configuration</b>	A-B-D-Ctrl-A1-H-A1-B / One profile per card
<b>Reagent</b>	Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping contains: <ul style="list-style-type: none"> <li>• Monoclonal Anti-A, Anti-B, Anti-AB and Anti-D (VI-) prefilled with gel from microtube 1-3</li> <li>• Microtube 4 contains neutral gel which serves as a negative control</li> <li>• Microtube 5 contains Anti-A1 lectin and Microtube 6 contains Anti-H lectin</li> <li>• Microtube 7 and 8 contains neutral gel for reverse grouping</li> </ul>
<b>Standardisation</b>	Follows FDA guidelines
<b>Storage / Stability</b>	4-25°C, 12 months

Presentation	Pack	Cat.No.
<b>Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping</b>	24 Cards	102560192



## Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping

**Intended Use:** \_\_\_\_\_

Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping is used for determination of ABO and Rho (D) blood group systems with subgrouping and reverse grouping for confirmation of ABO grouping.

### Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping kit components:

<b>Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping</b>	Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping contains Monoclonal Anti-A, Anti-B, Anti-D (VI-) prefilled with gel from microtube 1-3, Microtube 4 contains neutral gel which serves as a negative control, Microtube 5 contains Anti-A1 lectin, Microtube 6 contains Anti-H lectin and Microtube 7 & 8 contains neutral gel for reverse grouping
<b>Configuration</b>	A-B-D-Ctrl-A1-H-A1-B / One profile per card
<b>Other Accessories</b>	Card fitment and Package Insert

### Performance Characteristics:

**Standardization:** Follows FDA guidelines

FEATURES	BENEFITS
<b>8 Microtubes</b>	Forward grouping along with Anti-A1 and Anti-H Reverse grouping along with Autocontrol
<b>Stable End Reactions</b>	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
<b>Objective Interpretation of results</b>	Easy Interpretation of weak reactions Relatively skill free
<b>Small sample size</b>	More flexible More convenient
<b>Indigenously Manufactured</b>	Easy availability Long shelf life Another first from Tulip
<b>Tulip's Technical Support</b>	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards