

Technical Data Sheet Xylene histological grade

CE IVD *For In Vitro Diagnostic Use* *Manufacturer: Histo-Line Laboratories Srl*

Code Product	Packaging
3410/20 Xylene for histology and cytology	4 bottle of 5000 ml
3410/2500 Xylene for histology and cytology	4 bottle of 2500 ml
Code CND: W01030799	

Intended Use

Procedures of the pathology and cytology laboratory. Miscible with alcohols (ethanol, isopropanol, butanol), waxes and mounting means.

Prepared for installation of cyto-histological samples for optical microscopy.

General features

Primary container: bottle in polyethylene terephthalate (PET). Useful capacity of 5 liters or 2.5 liters. White color.

Polyethylene cap with leak-proof seal.

Secondary container: cardboard box.

Functional features

Components: Xylene mixture of isomers CAS 1330-20-7

Storage Condition

Keep the reagent at a lower temperature of 30 ° C. Keep containers tightly closed.

Validity

After first opening, the product is usable until the expiration date, if correctly stored. Validity: 3 years.

Disposal

Hazardous waste; confer on specialized and authorized companies, according to current legislation.

Carefully read the information on the classification of dangerous substances on the label. Always consult the MSDS where can I find information about the risks represented by the mixture, the precautionary measures during use, first aid measures and intervention in case of accidental release. Do not use if damaged primary container.

N.B. Although not reported in this technical sheet: compositions, precautionary measures, risk phrases, refer to storage, fire precautions, and more, refer to the Safety Data Sheet can be downloaded from our website: www.histoline.com or request it to: info@histoline.com

Milan Jul 2021