



	Neutral Gel Card
Application	For saline and Enzyme phase tests including antibody screening, antibody identification, compatibility testing and ABO reverse grouping
Card Configuration	6 x Neutral gel / Profile per card depends on test parameter
Reagent	Matrix™ ABO reverse grouping card contains: <ul style="list-style-type: none"> • Neutral gel in all six microtubes
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ Neutral Gel card	24 cards	102740024





Matrix™ Neutral Gel Card

Intended Use: _____

Matrix™ Neutral Gel Card is used for saline and enzyme phase tests including antibody screening, antibody identification, compatibility testing and ABO Reverse Grouping. Contains neutral gel.

Matrix™ Neutral Gel Card kit components:

Matrix™ Neutral Gel Card	Matrix™ Neutral Gel card contains six microtubes pre-filled with neutral gel in a suitable buffer and preservatives
Configuration	6 x neutral / Profiles per card depends on test parameter
Other Accessories	Card fitment and Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines

FEATURES	BENEFITS
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Enhanced Sensitivity	Better detectability especially in antibody screen and cross match
Eliminates cell wash step in IAT	Less skill and technique demands Saves time
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards