

# Human Immunodeficiency Virus Type 1 RNA Quantitative Fluorescence Diagnostic Kit

(PCR-Fluorescence Probing)

- **Dual gene targets**
- **Migh sensitivity**
- **V** Full automation enabled

## **EXTREMELY SIMPLE WORKFLOW**

1. Upload sample and reagent

2. Select program

3. Detect and analyze





# **FEATURES**









#### **HIGH SENSITIVITY**

25 IU/mL, meeting HIV infection clinic guidelines

### **GENOTYPE-WIDE** COVERAGE

HIV-1 group M, O, N

#### **ADVANCED EXTRACTION TECHNOLOGY**

Super-paramagnetic nano-beads, one wash, no heating

#### **FULLY AUTOMATED**

Matching Natch CS2, which automates nucleic acid extraction for standardised detection

#### **PARAMETERS**



Items	Parameter	
Specimen Type	Plasma	
Extraction Platform	Advanced magnetic beads technology	
Target Genes	HIV-1 RNA LTR and GAG gene	
Genotype	HIV-1 group M, O, N	
Internal Control	Exogenous IC	
PCR Instrument	Real-time PCR instruments with FAM,ROX and HEX/VIC channels(e.g. SLAN-96P, ABI7500, Roche cobas 480,etc.)	
Sensitivity	25 IU/mL	
Linear range	50 - 1.0E+08 IU/mL	

### **ORDER INFORMATION**

No.	Product Name	Spec.
S3046E	Human Immunodeficiency Virus Type 1 RNA Quantitative Fluorescence Diagnostic Kit (PCR-Fluorescence Probing)	48T/Kit
S1001E	Nucleic Acid (RNA) Extraction or Purification Kit (Magnetic beads method)	24T/Kit



ADDR.: No. 680, Lusong Road, Yuetu District, Changsha, Hunan Province, 410205, P.R.China

**TEL.**: +86-731-88883176-6116 **FAX**: +86-731-89736985

MAIL: info@sansure.com.cn / support@sansure.com.cn

WEB: www.sansureglobal.com

